


Methods for Wastewater Disease Surveillance: A Reagent Manufacturer's Perspective

Presented by:
Subhanjan Mondal, Ph.D.
Sr. Research Scientist



Core Technologies in:

- **Nucleic Acid Extraction and Amplification**
- **Protein Expression, Purification and Analysis**
- **Cell Health, Viability and Energy Metabolism**
- **Cellular, Genetic and Protein Reporters**



Our “Wastewater Surveillance Journey”

SARS-CoV-2
detected in
wastewater

Interlaboratory
Study
Participation ->

Direct Capture
Method
developed

Kits available in
catalogue

Automation Kits
including Ceres
Nanotrap® A

March 2020

April 2020

Spring 2020

July 2020

October
2020

Q2 2021

Nov 2021

Q1 2022

Ongoing

Adoption in
>50 countries

KWR

Global
Wastewater
Surveillance
Summit



RT-qPCR
Detection kit for
wastewater
including quant
standards and
controls




Method validation
paper in STOTEN



A direct capture method for purification and detection of viral nucleic acid enables epidemiological surveillance of SARS-CoV-2
Subhanjan Mondal ^{1,*}, Nathan Feirer ¹, Michael Brockman, Melanie A. Preston, Sarah J. Teter, Dongping Ma, Said A. Goueli, Sameer Moorji, Brigitta Saul, James J. Cali
Promega Corporation, 5400 E Cheryl Pkwy, Fitchburg, WI 53711, United States of America

Additional
methods and
qPCR Kits
under
development



Challenges of Wastewater Samples for Molecular Biology

- Sample composition greatly **varies** between locations/sewersheds
- Analytes are present in **low abundance**
- Wastewater contains compounds that **inhibit PCR**
- “**Solids**” vs “**Liquid**” phase – which one to use?
- Simple, scalable and consistent workflow is needed for surveillance.

Madison Metropolitan
Sewerage District



The Wastewater : A simplistic view



Dietary Fibers
(cellulose)



Tissue paper
(cellulose)




Particulates
(silica)

Chemicals


- Urea
- Detergents
- Salts

PCR inhibitors


- Bilirubin
- Humic acid
- Fulvic acid




Virus




Bacteria




Nucleic acid
naked



Nucleic acid
protected



The Inspiration: Sewage, Salt, Silica and SARS-CoV2 (4S) method



Sewage, Salt, Silica, and SARS-CoV-2 (4S): An Economical Kit-Free Method for Direct Capture of SARS-CoV-2 RNA from Wastewater

Oscar N. Whitney, Lauren C. Kennedy, Vinson B. Fan, Adrian Hinkle, Rose Kantor, Hannah Greenwald, Alexander Crits-Christoph, Basem Al-Shayeb, Mira Chaplin, Anna C. Maurer, Robert Tjian, and Kara L. Nelson*



Cite This: *Environ. Sci. Technol.* 2021, 55, 4880–4888



[Read Online](#)

Sample Concentration







Release TNA from solid fraction
30' Protease

Prepare TNA for capture on silica:
Chaotropic agents
+ alcohol

Bind TNA
to Midi-column
on Vacuum
Manifold
Wash & Inhibitor Removal

Elute TNA
1ml Eluate

0.5-1 ml TNA



Wizard® Enviro TNA Kit

Direct Capture & Spin column



Maxwell RSC® Enviro TNA Kit


Direct Capture & Maxwell® RSC Cartridges



Advantages of Direct TNA Capture

**Total
Nucleic Acid
concentration**


- 40 ml input provides higher yield than PEG/NaCl Method
- Consistent Recovery is ideal for surveillance



Recovery of **SARS-CoV-2** genetic material was more than 20-times more than the PEG/NaCl method.
Mondal et al, 2021

% Recovery of Direct Capture Method	
Spike virus*	% Recovery
MS2 phage	39.67 ± 10.66
OC43 coronavirus	63.13 ± 4.16
229E coronavirus	40.09 ± 10.36

- **Consistent yield**
- **Good Recovery**




High Yield and Quality for Sequencing

Case Report

Evaluation of viral concentration and extraction methods for SARS-CoV-2 recovery from wastewater using droplet digital and quantitative RT-PCR

Lampros Dimitrakopoulos¹, Aikaterini Kontou¹, Areti Strati, Aikaterini Galani, Marios Kostakis, Vasileios Kapes, Evrikleia Lianidou, Nikolaos Thomaidis^{**}, Athina Markou^{*}

Laboratory of Analytical Chemistry, Department of Chemistry, National and Kapodistrian University of Athens, University Campus, Zografou, 15771, Athens, Greece



Recoveries (%) of the SARS-CoV-2 spiked viral standard (standard deviation values (SD)) on three different days with all five different concentration methods and all four different extraction kits based on RT-qPCR.

Concentration method	DAY 1 Recovery (SD)	DAY 2 Recovery (SD)	DAY 3 Recovery (SD)	AVERAGE
PEG precipitation	0.1289 (0.0165)	0.1067 (0.0123)	0.0864 (0.0142)	0.1073
PEG & Glycine	0.1373 (0.0082)	0.2703 (0.0217)	0.3031 (0.0042)	0.2369
Ultrafiltration	0.0104 (—)	0.0018 (0.0000)	0.0115 (0.0065)	0.0079
Adsorption-extraction	0.0202 (0.0044)	0.0042 (0.0028)	N/A	0.0122
Direct capture	1.82 (0.06)	1.63 (0.01)	2.73 (0.02)	2.06
Extraction method	DAY 1 Recovery (SD)	DAY 2 Recovery (SD)	DAY 3 Recovery (SD)	AVERAGE
IDEXX	6.31 (0.0012)	7.07 (0.0039)	12.7 (0.0054)	8.7
QIAGEN Viral	22.0 (0.0035)	14.8 (0.013)	25.6 (0.0122)	20.8
QIAGEN Water	1.02 (0.0016)	0.90 (0.0013)	4.41 (0.0028)	2.11
Promega	39.3 (0.0071)	36.0 (0.0239)	37.6 (0.0258)	37.6

Enviro TNA kits yield more nucleic acid from a smaller sample volume than other methods

Evaluation of two different concentration methods for surveillance of human viruses in sewage and their effects on SARS-CoV-2 sequencing




Inés Girón-Guzmán^{a,1}, Azahara Díaz-Reolid^{a,1}, Enric Cuevas-Ferrando^a, Irene Falcó^a, Pablo Cano-Jiménez^{b,c}, Iñaki Comas^{b,c}, Alba Pérez-Cataluña^{a,*}, Gloria Sánchez^a

^a Department of Preservation and Food Safety Technologies, Institute of Agrochemistry and Food Technology, IATA-CSIC, Av. Agustín Escardino 7, Paterna 46980, Valencia, Spain

^b Instituto de Biomedicina de Valencia (IBV-CSIC), C/ Jaime Roig, 11, Valencia 46010, Spain

^c CIBER in Epidemiology and Public Health (CIBERESP), Valencia, Spain




TNA purified using AICl3 ppt

TNA purified using Maxwell® Enviro TNA

Enviro TNA purified nucleic acid provides high quality sequencing data

New Purification Methods: Going after the “solids”



APPLIED AND ENVIRONMENTAL MICROBIOLOGY, Mar. 1976, p. 354-358
Copyright © 1976 American Society for Microbiology

Vol. 31, No. 3
Printed in U.S.A.

Demonstration of Solids-Associated Virus in Wastewater and Sludge

FLORA MAE WELLINGS,* ARTHUR L. LEWIS, AND CARROL W. MOUNTAIN
Epidemiology Research Center, State of Florida, Department of Health and Rehabilitative Services, Tampa, Florida 33614


Comparative analysis of Adsorption-Extraction (AE) and Nanotrap®
Magnetic Virus Particles (NMVP) workflows for the recovery of endogenous
enveloped and non-enveloped viruses in wastewater

Warish Ahmed ^{a,*}, Aaron Bivins ^b, Asja Korajkic ^c, Suzanne Metcalfe ^a, Wendy J.M. Smith ^a, Stuart L. Simpson ^d


Adsorption of Respiratory Syncytial Virus, Rhinovirus, SARS-CoV-2, and F+ Bacteriophage MS2 RNA onto Wastewater Solids from Raw Wastewater

Laura Roldan-Hernandez and Alexandria B. Boehm*

 Cite This: *Environ. Sci. Technol.* 2023, 57, 13346–13355

 Read Online

Detection and Quantification




GoTaq® Enviro qPCR kits
made for **wastewater surveillance**:

- PCR inhibitor tolerant qPCR and RT-qPCR Master-Mix
- Multiplexed to include process control (PMMoV or CrAssphage)
- **Quantitation standard RNA/DNA** included



Design of the Wastewater SARS-CoV-2 RT-qPCR Detection Kit




Detection Channel in Multiplex				
Cat.#	FAM	HEX	ROX	Cy5
AM2150	N1	N2	IAC	PMMoV
AM2160	E	N2	IAC	PMMoV
AM2100-2130	N1, N2, or E	IAC		PMMoV
CS317431	Flu A	Flu B	SC2	PMMoV

Quantitation Standard RNAs are available for: SARS-CoV-2 (N+E), SC2, PMMoV, FluA, FluB


Assays with Proper QA/QC controls

Quantitation Standards



Quantify RNA using QuantiFluor RNA System and droplet digital RT-PCR.


Dilute the RNA at 4×10^6 copies/ul in TE-Buffer + protectant



Human Fecal Control

Pepper mild mottle virus

- An **RNA** plant virus that is commonly found in human feces (diet)
- It is present in wastewater globally and serves as **human fecal indicator**



PMMoV levels in wastewater from one sampling point over multiple days provided by Zeynep Cetecioglu, KTH, Stockholm


Internal Amplification Control

- An exogenous RNA that is in every qPCR assay serving as an amplification control
- Template for the Internal control used is Firefly Luciferase RNA

A shift in Cq value from the no template control (NTC) reflects the level of RT-qPCR inhibition in a sample.

$$\Delta Cq = Cq_{[Sample]} - Cq_{[NTC]}$$

$\Delta Cq > 2$ represents significant inhibition of the reaction




Application scientists develop additional applications:

[illegible]

Promega purification
used for qPCR and for
**Oxford Nanopore
Sequencing**

[illegible]

Automation for Concentration and purification




Product Application

Automated Total Nucleic Acid Purification from 40ml of Wastewater using Streamlined Protocol

Purify total nucleic acid (TNA) from 40ml of wastewater with the Maxwell® RSC Enviro TNA Kit using a streamlined protocol to reduce preparation and processing steps.

<p>KIT: Maxwell® RSC Enviro TNA Kit (Cat #A51831)</p> <p>Analyses: RT-qPCR</p> <p>Sample Type(s): Wastewater</p> <p>Input: 40ml</p> <p>Materials Required:</p> <ul style="list-style-type: none"> * Maxwell® RSC Enviro TNA Kit (Cat #A51831) * Isopropanol * 95% Ethanol * 250ml conical tubes * TubeRite centrifuge capable of spinning 250ml conical tubes, 3000 x g * Vacuum Pump * Vac-Max™ Laboratory Vacuum Manifold (Cat #A7231) * Elututor™ Vacuum Ejection Device (Maddox) * Maxwell® RSC Instrument (Cat #AS4500 or Maxwell® RSC 48 Instrument (Cat #AS8500) * Heat block or waterbath set to 60°C 	<div style="border: 1px solid black; padding: 10px;"> <p>This protocol was developed by Promega Applications Scientists and is intended for reference use only. Users are responsible for determining suitability of the protocol for their application. For further information, see Technical Manual TM0003, available at www.promega.com/products-and-services/technical-manuals/ or contact Technical Support at support@promega.com.</p> </div>
--	---



Protocol:

1. Prepare Cocktail Wash 1 & 2 as indicated on the bottles.
2. Combine 12ml of Binding Buffer 1 and 1 ml of Binding Buffer 2 for each sample* Cocktail A (BCA).

*Alternatively, a master mix may be prepared, up to least 200x coverage mix of BCA with 20% overage; prepare additional Binding Buffer 1 & 2 purchased. The volumes provided in the Maxwell® RSC Enviro TNA kit at 47 reactions when considering the preparation of 50% overage.

3. Preheat 600µl of Nuclease-Free Water, per sample, to 60°C.
4. Capture and Concentration:
 - a) Collect 40ml of wastewater in a 250ml conical tube.
 - b) Add 5.5mL of Proteinase Solution. Insert mix and incubate for 30 minutes at 47°C.
 - c) Add 13ml of BCA.

PROMEGA CORPORATION • 8000 REDWOOD BLVD • MADISON, WI 53704 USA • TEL: 608.276.8000 FAX: 608.276.8001
 SCIENTIFIC APPLICATIONS NOTES PAGE 1110 - 03/2017

Reduce hands-on time by using larger containers

[illegible]

Compare Standard RNAs

[illegible]

Using additional qPCR instruments


[illegible]

...and more

Where Next?

Now that the COVID-19 pandemic is behind us, how do we sustain the surveillance system that was set up during the pandemic

- Targets?
- What sample-type to use that's best for the targets
- What throughput



Please stop by our booth
Reach us at:
applied@promega.com



Thank you